


CETSA®EXPLORE: Target Engagement Profiles to Study SAR, Off-Targets, and Toxicity at Proteome Wide Level



IM CABALLERO ; T FRIMAN; A CHERNOBROVKIN; E GILSON; S LUNDGREN; D MARTINEZ MOLINA; PELAGO BIOSCIENCE AB, SWEDEN




TARGET ENGAGEMENT ATLAS OR TEA

The objective of this project is to collect compound target engagement profiling, or CETSA® fingerprints for a wide spectrum of compounds representing different parts of the chemical, pharmacological and biological space.

Positive unintended effects, as well as liabilities can be discovered at an early phase of drug development.


TARGET ENGAGEMENT ATLAS PIPELINE



TARGET ENGAGEMENT PROFILES FOR SAR STUDIES

Some chemical motifs have a strong impact on the CETSA fingerprint, which can be applied to SAR studies.

Carbonic Anhydrase 13 profiling shows its inhibitor Acetazolamide as a hit, as well as some coxibs, which are NSAID class, and MoA reported to act through inhibition of cyclooxygenase 2. All these compounds had a sulfonamide group in common, and this chemical feature was not found in a Coxib which didn't shift Carbonic anhydrase 13.




TEA AND THE CELLULAR CONTEXT

Target engagement atlas identifies hallmarks of cellular responses to different types of drugs and assigns biomarkers beyond specific target engagement.

We chose three different cell lines that cover most of the human proteome: HepG2, (human hepatocellular carcinoma), U87MG (primary glioblastoma), and K562 (myeloid leukemia, lymphoblast).

A large number of compounds proteomic profiles correlate with each other, primary target in particular, but other targets are cell line-specific. We observe clear cellular discrepancies in response patterns to a large number of drugs




Visualization of the results with tSNE in HepG2 cells, an unsupervised non-linear reduction model that group related compound fingerprints close to each other. Colour coded by pharmacology, compounds tend to cluster together within their pharmacological effect.


CONCLUSIONS

CETSA® EXPLORE compound profiling of hundred of known drugs and probes in three different cell lines reveals:

- Specific chemical motifs influence compound CETSA® fingerprints.
- Compound protein hits response differently dependent on the cellular background.



Scatter plot to correlate proteomic responses in K562 and HepG2 or K562 and U87MG with 100 compounds .



Different Src type kinases response differently to treatment in different cell types.

Yes1 and Src responses correlated between HepG2 and U87MG but not in K562.

Dasatinib, a specific SRC inhibitor, is effective in the three cell lines. Interestingly, Dasatinib is used in clinics for bcr-abl Leukemia, the same source tumor as K562.

- References:
- Martinez Molina Science 2013
 - Savitski Science 2014
 - Chernobrovkin BioRx 2020
 - Friman Bioorg. Med. Chem 2020
 - Friman SLAS Discovery 2021
 - Chernobrovkin SLAS Discovery 2021
 - Hendricks ACS Chemical Biology 2022