

In Vitro Drug Effects on Cancer Cell Morphology and Functional State Revealed by Multiparameter Imaging Mass Cytometry

Olga Ornatsky, Alexandre Bouzekri, Roberto Spada, Bedilu Allo and Jessica Watson

FLUIDIGM®

Introduction

Results of *in vitro* cell culture drug testing are correlated with clinical response to chemotherapy: The benefit of using *in vitro* models lies in the ability to probe cellular response in a controlled closed system, where the effects of drug concentrations, treatment duration, drug efflux kinetics and multidrug combinations can be assessed using a variety of cell biology techniques. However, obtaining highly multiplexed data for surface, intracellular and functional markers for each cell within a test well has been challenging using standard immunofluorescence techniques. Imaging mass cytometry (IMC) (1,2) allows us to investigate complex effects of chemotherapy drugs as well as the intracellular localization of metal-containing drugs (e.g., cisplatin and oxaliplatin), ruthenium-containing drugs (e.g., NAMI) and cytostatic agents (e.g., nocodazole and etoposide) (Table 2). A large panel of metal-tagged antibodies (Table 1) was used to analyze responses to drug treatments at the single-cell level. Proteins involved in DNA damage repair, apoptosis, cell proliferation, substrate adhesion, organelle morphology and signaling pathways, as well as surface and cytoskeletal markers, were studied. S-phase cells were visualized by detection of iodine (¹²⁷I) from 5-iodo-2'-deoxyuridine (IdU) incorporated by growing cells.

Imaging mass cytometry workflow. The immunostained and dried samples of tissue sections (FFPE or cryosections) and cells attached to glass slides are inserted into the ablation chamber of the Hyperion™ Imaging System (Fluidigm), where a 1 µm spot-size pulsed UV laser ablates the tissue. Isotopes associated with each spot are transferred into the mass cytometer, detected and indexed against the source location, yielding an intensity map of the target proteins throughout the tissue.




Table 1. Metal-conjugated antibodies against human antigens


Metal-conjugated antibodies to surface markers	Clone	Fluidigm/other vendor catalog number	Metal-conjugated antibodies to intracellular markers	Clone	Fluidigm/other vendor catalog number
B7-H3- ¹²⁵ Yb	D9M2L	CST™ 140588F	Beta-Catenin- ¹⁴⁴ Sm	D10A8	3147005A
CD44- ¹⁴⁴ Nd	NY2	BioLegend® 365202	Beta-Tubulin- ¹²⁷ Yb	rabbit poly	Abcam ab6046
CD24- ¹⁵⁶ Er	ML5	3169004B	CD107a- ¹⁵² Eu	H4A3	3151002B
CD326- ¹⁴⁶ Pr	9C4	3141006B	c-Myc- ¹⁵⁷ Er	9E10	BioLegend® 626802
CD47- ²⁰³ Bi	CC2C6	3209004B	cPARP(Asp214)- ¹⁴⁵ Nd	F21-852	3143011A
CD44- ¹⁶⁶ Er	BJ18	3166001B	CyclinB1- ¹⁶² Dy	GNS-1	3164010A
CD49e- ¹⁶² Gd	H110a	3156001B	Histone3- ¹⁵⁵ Yb	DIH2	3176016A
CD51- ¹⁶⁵ Tm	NY2	BioLegend 327902	Keratin (CK8/18)- ¹²⁷ Yb	C51	3174014A
CD59- ¹⁷² Tb	p282 (H19)	3173009B	Ki67- ¹⁶⁴ Er	B56	3168007B
CD61- ¹⁵² Nd	VI-PL2	3150001B	p53- ¹⁵² Nd	DO-7	Abcam ab80644
CD71- ¹⁷⁵ Lu	OKT-9	3175011B	Pan-Keratin- ¹⁶² Dy	C11	3162027A
CD9- ¹⁷² Yb	SN4 C3-3A2	3171009B	pH2AX(pS139)- ¹⁶⁵ Ho	N1-431	BD Biosciences 560443
CD98- ¹⁶⁵ Tb	H4A3	3151002B	pHistone3- ¹⁷² Lu	Ser28 (HTA28)	3175012A
ErbB2/HER2- ¹⁵² Eu	24D2	BioLegend 324402	Thioredoxin- ¹⁴⁴ Nd	2G11/TRX	3146016B
EGFR- ¹⁷² Er	AY13	3170009B	Vimentin- ¹⁵⁸ Gd	RV202	3156023A
			Intercalator-103Rh	N/A	201103B
			Intercalator-Ir	N/A	201192B
			IdU	N/A	201127

Table 2. Drugs and cell lines



Drug	Amount	Drug treatment	Drug effect	Cell line tested
Cisplatin	5 µM	24 hours	DNA damage	A-431 (ATCC® CRL-1555™) Skin epidermoid carcinoma
Cisplatin	20 µM	24 hours	DNA damage	HeLa (ATCC CCL-2™) Cervix adenocarcinoma
NAMI	25 µM	2 hours	Migration and adhesion inhibitor, cell cycle arrest, induction of apoptosis	A-431 (ATCC CRL-1555) Skin epidermoid carcinoma
Etoposide	10 µM	18 hours	DNA damage	HeLa (ATCC CCL-2) Cervix adenocarcinoma
Nocodazole	10 ng/mL	24 hours	Cell cycle arrest in mitosis	HeLa (ATCC CCL-2) Cervix adenocarcinoma

Results

IMC false-color images and quantitative analysis of SKOV3



Cells were grown on chamber slides (Lab-Tek™ II Chamber Slide™) for several days and stained with a mixture of metal-conjugated antibodies. SKOV3 cell culture 16-bit TIFF images were generated from raw IMC files (MCD files). Scale bar = 100 µm. Nuclei were stained with DNA-binding Ir-intercalator. The Ir-intercalator DNA ¹⁹¹/¹⁹³Ir dot plot shows correlation of the intensity of the two channels as a red density of points corresponding to single cells.




SKOV3 cell heterogeneity analyzed by viSNE: (A) nuclear markers, (B) surface markers, (C) cytoplasmic markers


4,229 segmentation areas corresponding to single cells were converted into an FCS file with mean values per object for each channel. 4,000 single cells were mapped to a two-dimension viSNE map, level of expression of each marker is shown using a heat map gradient. Analysis of the FCS file generated from raw data (not normalized) is then performed using the Cytobank platform. Using single-parameter histograms, comparative dot plots and viSNE analysis (4), we can map the multiparametric data from an IMC region into different formats.

In vitro drug testing on adherent cell lines


A431 cell culture



NAMI, 25 µM 2 hr



HeLa cell culture



Etoposide is a topoisomerase inhibitor. Etoposide forms a ternary complex with DNA and topoisomerase II enzyme, prevents re-ligation of the DNA strands and by doing so causes DNA strands to break. Cancer cells rely on this enzyme more than healthy cells, since they divide more rapidly. This causes errors in DNA synthesis and promotes apoptosis of cancer cells. Etoposide is used as a form of chemotherapy for cancers such as sarcomas, lung, testicular cancers, lymphoma and nonlymphocytic leukemia and glioblastoma.

Nocodazole is an antineoplastic agent that exerts its effect in cells by interfering with the polymerization of microtubules. Nocodazole stimulates the expression of LATS2, which potentially inhibits the Wnt signaling pathway. Microscopy of nocodazole-treated cells shows that they do enter mitosis but cannot form metaphase spindles because microtubules (of which the spindles are made) cannot polymerize.

References

- Chang, Q., Ornatsky, O., Siddiqui, I. et al. "Biodistribution of cisplatin revealed by imaging mass cytometry identifies extensive collagen binding in tumor and normal tissues". *Scientific Reports* 6 (2016): 36641.
- Chang, Q., Ornatsky, O., Siddiqui, I. et al. "Imaging mass cytometry." *Cytometry A* 91 (2017): 160–169.
- de Chaumont, F., Dallongeville, S., Chenouard, N. et al. "Icy: an open bioimage informatics platform for extended reproducible research." *Nature Methods* 9 (2012): 690–696.
- Amir el-AD, Davis, K.L., Tadmor, M.D. et al. "viSNE enables visualization of high dimensional single-cell data and reveals phenotypic heterogeneity of leukemia." *Nature Biotechnology* 31 (2013): 545–552.