

An effort to build an in vitro high-throughput screening for seizure liability based on calcium oscillation of human iPSC-derived dopaminergic neurons with astrocyte

Sunao Hisada¹, Yiwen Wang², Du Shouming², Natsumi Kato¹, Masanori Matsubara¹ and Fumio Iwase¹ Systems Division, Hamamatsu Photonics K.K., ¹Hamamatsu, Japan, ²Bridgewater, NJ, USA
Aoi Odawara³ and Ikuro Suzuki³ ³Department of Electronics and Intelligent Systems, Tohoku Institute of Technology, Miyagi, Japan


1. Introduction

Drug-induced seizure liabilities are typically tested in expensive and low throughput ex vivo rat hippocampal brain slices. An in vitro model based on micro-electrode arrays using human iPSC derived neurons has been suggested (<https://doi.org/10.1093/toxsci/kfy029>). There is a surge of interest in higher throughput methods, one of them being fluorescence measurement of calcium oscillations in neuronal cultures. However, progress has been limited by the long imaging acquisition time for an entire 96- or 384-well plate, and also by the heterogeneity of neurons derived from human iPSCs. With the advent of new camera technologies as well as availability of highly enriched, functionally mature human neurons derived from iPSCs, these limitations can now be overcome. Here we attempt to build a model to compare human iPSC derived neurons and astrocyte, with rat primary cells upon treatment of reference and control compounds.


Here, we optimized the assay conditions for a calcium-sensitive fluorescent dye (Cal-520AM) to Ca²⁺ transients in respective neurons. The fluorescence measurements for a 384 plate were performed with 70 ms interval for 5 min at 37°C, in the high-throughput manner, using a camera-based kinetic plate reader, FDSS (Hamamatsu), which measures fluorescence signals of all 384 wells in a microplate simultaneously. The burst rates of neurons (Primary Rat neuron and Cellular Dynamics International iCell® Dopa neuron and iPSC astrocyte) labeled with Cal-520AM, which were estimated from fluorescence signal. Then, we examined effects of 11 reference compounds including 4-AP, Picrotoxin on Ca²⁺ transient-waveforms to detect drug-induced waveform changes and burst rate changes, etc. Such high-throughput fluorescence measurements of Ca²⁺ transient in primary or hiPSC derived neurons would be greatly useful for assessment of the seizure risk of compounds, particularly at the early stages of drug discovery..



FDSS/uCELL



200Hz Ca²⁺ oscillation data



dispenser

384/96-well plate

EM-CCD camera

2. Material and Method

Compound	Function	Dose1	Dose2	Dose3	Dose4	Dose5
4-Aminopyridine(4-AP)	Seizurogenic Compounds – Potassium Channel Blocker	>0.3	1	>3	10	>30
Pentylenetetrazole(PTZ)	GABA antagonist,circulatory and respiratory	10	>30	100	>300	1000
Picrotoxin (PTX)	Toxin(GABA antagonist-Seizurogenic Compounds)	0.1	>0.3	1	>3	10
Linopirdine	acetylcholine releaser(CNS);release of acetylcholine and other neurotransmitters	1	>3	10	>30	100
Pilocarpine	eye drop use(cholinergic)	>0.3	1	>3	10	>30
Strychnine	Glycine Receptor Antagonist-Seizurogenic Compounds)	>0.3	1	>3	10	>30
Amoxapine	antidepressant(CNS);inhibits 5-HT2 serotonergic receptors	>0.3	1	>3	10	>30
Chlorpromazine(CPZ)	antipsychotic(CNS)	0.1	>0.3	1	>3	10
Phenytoin	anti-seizure	1	>3	10	>30	100
Amoxicillin	antibiotic	1	>3	10	>30	100
Acetaminophen	Negative compounds	1	>3	10	>30	100
DMSO	Vehicle	0.1	(-)	(-)	(-)	(-)

Sample preparations


Rat primary neuron(E15), iCell® Dopa neuron and iPS Astrocytes (Cellular Dynamics International) were plated to the PEI-coated 384-wellplates (Corning) or Nanofiber plate(Nanofiber Solutions (Columbus, OH)). at 35,000 cells/well and incubated in 80 ul of the culture medium. On the day 14-21 after plating, the cells were loaded with fluorescent dyes for Ca²⁺ transient measurements. For measurements of Ca²⁺ transients, 50 µg of the dye (Cal-520 AM (AAT Bioquest)) was diluted in 45 µl DMSO and then in 20 ml of culture medium (CDI) containing 1.25 mM probenecid (Sigma-Aldrich) warmed at 37 °C (2 µM of dye final concentration). After the culture medium was removed, 72 ul of the dye solution was applied to the well, and incubated for 45 min at 37 °C in 5% CO₂. The measurements were performed with low Mg²⁺ buffer(0.1mM Mg²⁺). .

After one fluorescence measurement in the absence of test compounds (pre), 8 ul of each test compound (10 x final concentration) in PBS (-) was applied to the wells, and the cells were incubated for 10-45 min at 37 °C in 5% CO₂. In this study, a library of pharmacologically active compounds (Office for University-Industry Collaboration, Osaka University, Japan) was used.


Fluorescence measurements and analysis

The fluorescence intensity changes of all wells in a microplate were taken simultaneously with every 70 ms for 5 min (4000 frames) at 37 °C using a kinetic plate reader, FDSS 7000EX or FDSS/uCELL (Hamamatsu). The intracellular Ca²⁺ concentration changes (calcium waveforms) were analyzed using the FDSS Waveform Analysis Software (Hamamatsu), which estimates peak rate, peak-to-peak time, peak width at 10%-90%, rising slope, falling slope, and more. All compounds were tested in N=4~5, 5 dose conditions. The result of each well was normalized with that of control (0.1 % DMSO).

3. FDSS Ca2+ oscillation data




FDSS Ca²⁺ oscillation image of Rat cortex neuron from 384 well Nanofiber plate




Ca²⁺ oscillation image of human iPSC-derived Dopa/Astrocytes co-cultured neuron (CDI) from 384well plate

4. Result


Spotfire analysis of FDSS drug-induced synchronization of Ca²⁺ oscillation data in Rat cortex neuron and human iPSC-derived Dopamine/Astrocytes co-cultured neurons (N=4~5, 5 dose, 384 well plate and using low Mg²⁺ buffer(0.1mM Mg²⁺))



Burst per 5 minutes (% DMSO Ctrl)



RFU Amplitude (%DMSO Ctrl)



RFU Area under curve (% DMSO Ctrl)

5. Result

Summary of effects of some reference compounds on Ca²⁺ transient-waveforms measured with fluorescent dyes in Rat Primary cortex neurons (E15) and human iPSC-derived neurons.

Compounds	Human clinical dose	Rat cortical neuron	Rat primary cortex neuron(0.1mM Mg ²⁺)			CDI iPS Dopa/astrocytes co-cultured (0.1mM Mg ²⁺)			
			BPM	AMP	AUC	Comments	BPM	AMP	AUC
4-Aminopyridine(4-AP)	seizure	burst coordination improved	↑ ↑ ↑			BPM over 25% up at 1uM	↑ ↑ ↑	↑ ↑ ↑	↑ ↑ ↑
Pentylenetetrazole(PTZ)	seizure	activity increased/Excitatory	↑ ↑			BPM over 10% up at 100uM			
Picrotoxin (PTX)	seizure	burst coordination improved/Excitatory	↑ ↑ ↑	↑	↑	BPM over 15% up at 0.3uM	↓		
Linopirdine	NA	NA	↑ ↑ ↑			BPM over 20% up at 3uM	↑ ↑ ↑	↓	
Pilocarpine	seizure	no effect or activity slightly increased/Excitatory	↑	↑	↑	AUC over 20% up at 10uM	↑ ↑ ↑	↓	
Strychnine	seizure	no effect or activity reduced	↑ ↑ ↓		↓	BPM over 15% up at 1uM	↓	↓ ↓ ↓	↓ ↓
Amoxapine	sedation	no effect or activity reduced	↑ ↑ ↓			BPM over 20% up at 1uM	↑ ↑ ↓	↓ ↓	↓ ↓
Chlorpromazine(CPZ)	sedation	activity reduced/inhibitory	↑ ↑ ↓	↓		BPM >-10% Inhibitory (1uM)	↓ ↓ ↓	↓ ↓ ↓	↓ ↓
Phenytoin	inhibitory	activity reduced	↓	↓		AMP >-10% Inhibitory (1uM)	↓ ↓ ↓	↓ ↓ ↓	↓ ↓
Amoxicillin	NA	NA	↑	NC	NC	BPM over 10% up at 30uM	NC	NC	NC
Acetaminophen	NA	NA	↑	NC	NC	no change	NC	NC	NC
DMSO	NA	NA	NC	NC	NC	no change	NC	NC	NC

6. Acknowledgement

We thank Dr Yoon-Jeong Kim, Dr Bang Zhong Lin, Dr Kazuto Nunomura, Dr Tsutomu Sasaki (Osaka university) , Dr. Ko Zushida (Cellular Dynamics International Japan) for their great support and helpful discussion on the experiments using human iPSC-derived Neuron. For further information, questions and comments, please contact shisada@sys.hpk.co.jp or SDU@hamamatsu.com www.hamamatsu.com/us/en/community/fdss/index.html **Reference** PLoS One. 2014 Jan 8;9(1):e84755